

CardioDetect User Manual

Comprehensive Guide for Patients, Doctors, and Administrators

Version: 1.0

Date: December 2025

Table of Contents

1. Introduction
 2. Getting Started
 - Registration & Login
 - Password Recovery
 3. Patient User Guide
 - Dashboard Overview
 - Running a Risk Assessment
 - Viewing History
 - Managing Your Profile
 4. Doctor User Guide
 - Doctor Dashboard
 - Patient Management
 - Uploading Medical Reports (OCR)
 - Analyzing Results
 5. Administrator User Guide
 - Admin Panel Overview
 - User Management
 - Approving Profile Changes
 - System Analytics
 6. Troubleshooting & Support
-

1. Introduction

Welcome to **CardioDetect**, an advanced AI-powered platform designed to assess cardiovascular disease risk early and accurately. This manual will guide you through the features specific to your role.

Key Features:

- * **AI Risk Prediction:** Instant analysis using advanced machine learning.
- * **Medical Report OCR:** Upload PDF/Image reports for automatic data extraction.
- * **Secure Platform:** HIPAA-compliant security with role-based access.

2. Getting Started

Registration & Login

1. Navigate to the **Home Page**.
2. Click **“Get Started”** or **“Register”**.
3. Fill in your details:
 - **Name:** Full legal name.
 - **Email:** Valid email address (used for verification).
 - **Password:** Strong password (min 8 chars, mixed case, numbers).
 - **Role:** Select “Patient” or “Doctor” (Doctor accounts require admin verification).
4. Click **“Create Account”**.
5. Check your email for a **Verification Link** and click it to activate your account.

Password Recovery

If you forget your password: 1. Click **“Forgot Password?”** on the Login screen. 2. Enter your registered email. 3. Click the link in the password reset email. 4. Set a new secure password.

3. Patient User Guide

As a patient, you can monitor your heart health, run assessments, and track changes over time.

Dashboard Overview

Your command center showing: * **Latest Risk Score:** Your most recent heart health assessment. * **Quick Actions:** Buttons to “New Assessment” or “View History”. * **Health Trends:** Simple charts showing your progress.

Running a Risk Assessment

1. Click **“Assess My Risk”**.
2. **Manual Entry:** Fill in the health form (Age, Blood Pressure, Cholesterol, etc.).
 - *Tip:* If you don’t know a value, use the default average provided or leave it blank if allowed.
3. Click **“Calculate Risk”**.
4. **View Results:**
 - **Risk Level:** Low (Green), Moderate (Yellow), High (Red).
 - **Probability:** Percentage chance of heart issues in the next 10 years.

- **Recommendations:** Personalized lifestyle advice based on your inputs.
5. **Download Report:** Click “Download PDF” to save a copy for your records.

Viewing History

- Click “**History**” to see a list of all your past assessments.
- Use the **Trends Graph** to see if your risk is increasing or decreasing.
- Click on any past entry to view full details.

Managing Your Profile

- Go to “**Profile**” to update personal details (Phone, Address).
 - **Sensitive Changes:** Changing critical health data (like Date of Birth) may require **Admin Approval** for security. You will be notified when approved.
 - **Data Rights:** You can request a “Data Export” or “Account Deletion” from the Settings tab (GDPR compliance).
-

4. Doctor User Guide

Doctors have powerful tools to manage patients and analyze medical documents efficiently.

Doctor Dashboard

- **Patient List:** Quick access to your assigned patients.
- **Recent Uploads:** Status of processed medical reports.
- **Risk Alerts:** Notifications for patients flagged as “High Risk”.

Patient Management

1. Go to “**My Patients**”.
2. **Search:** Find patients by name or ID.
3. **Assign Patient:** Use the “Add Patient” button to link a patient to your care (requires patient consent/code).
4. **View Profile:** Click a patient to see their full assessment history and trends.

Uploading Medical Reports (OCR)

Eliminate manual data entry by uploading scanned reports.

1. Click “**Upload Report**”.

2. Select a patient.
3. **Drag & Drop** a file (PDF, JPG, PNG).
 - *Supported:* Lab reports, clinical summaries.
 - *Max Size:* 10MB.
4. Click **“Analyze”**.
5. **Review Extracted Data:** The AI will extract values (Cholesterol, BP, etc.).
 - *Verify:* Check the “Confidence Score”.
 - *Edit:* Correct any misread numbers before submitting.
6. **Confirm:** The system runs the risk prediction automatically using the extracted data.

Analyzing Results

- **Clinical Report:** Generated PDFs include:
 - Patient Demographics.
 - Extracted Vitals.
 - AI Risk Assessment.
 - **SHAP Analysis:** “Why this prediction?” (Explains top contributing factors like ‘High Systolic BP’).
-

5. Administrator User Guide

Admins oversee the platform’s security, user base, and data integrity.

Admin Panel Overview

Access the comprehensive dashboard at `/admin`. * **System Stats:** Total users, daily assessments, error rates. * **Pending Actions:** Approval requests, user verifications.

User Management

- **View Users:** List all registered Patients and Doctors.
- **Actions:**
 - **Deactivate:** Suspend suspicious accounts.
 - **Reset MFA:** Assist users locked out of accounts.
 - **Assign Role:** Upgrade a user to “Doctor” or “Admin” status.

Approving Profile Changes

For data integrity, sensitive profile updates (e.g., Name, DOB) enter a “Pending” state. 1. Go to **“Approvals”**. 2. Review the **“Pending Changes”** list. 3.

Compare “Old Value” vs. “New Value”. 4. Click “**Approve**” to update the potentially or “**Reject**” to discard.

System Analytics

- **Prediction Logs:** detailed audit trail of every prediction run.
- **Performance Metrics:** Monitor API latency and ML model accuracy.
- **Security Logs:** View failed login attempts and potential brute-force attacks.

6. Troubleshooting & Support

Issue	Possible Cause	Solution
Login Failed	Incorrect password or unverified email	Check “Caps Lock” or check email for verification link.
OCR Reading Error	Low quality image or handwriting	Use a scanned PDF or clear photo. Ensure text is typed, not handwritten.
“Access Denied”	Insufficient permissions	Contact Admin if you need Doctor access.
Report Download Fails	Pop-up blocker	Allow pop-ups for CardioDetect.

Still need help? Contact support at support@cardiodetect.com.