

# CardioDetect Milestone 2 Report

---

## AI-Powered Cardiovascular Disease Risk Prediction System

---

**Project:** Early Detection of Heart Disease Risk

**Version:** 2.0

**Date:** December 2025

---

## Table of Contents

---

1. [Executive Summary](#)
  2. [Data Quality & Preprocessing](#)
  3. [Model Architecture](#)
  4. [Classification Models Comparison](#)
  5. [Regression Models Comparison](#)
  6. [Hyperparameter Tuning](#)
  7. [Model Evaluation Metrics](#)
  8. [Risk Categorization System](#)
  9. [Clinical Override Rules](#)
  10. [Saved Models & Pipelines](#)
  11. [Testing & Validation](#)
  12. [Conclusion](#)
- 

## 1. Executive Summary

---

CardioDetect is an AI-powered system for early detection of cardiovascular disease risk. This milestone delivers:

- **18 Classification Models** trained and compared

- **4 Regression Models** for continuous risk prediction
- **Production-ready MLP Classifier** achieving **99.25% accuracy**
- **Clinical Override System** to catch edge cases missed by ML
- **End-to-end OCR Pipeline** for medical document processing

Key Achievements

Metric	Target	Achieved
Classification Accuracy	>85%	99.25%
Model Precision	High	99.6% (MLP Binary)
Model Recall	High	99.6% (MLP Binary)
ROC-AUC Score	>0.95	0.999 (MLP Binary)

2. Data Quality & Preprocessing

Dataset Overview

Split	Samples	Percentage
Training	11,286	70%
Validation	2,418	15%
Test	2,419	15%
Total	16,123	100%

Data Source


- **Primary:** Framingham Heart Study dataset
- **Features:** 14 clinical measurements + engineered features
- **Target:** 3-class cardiovascular risk (LOW/MODERATE/HIGH)

## Missing Value Handling

Feature	Missing %	Imputation Method
guideline_risk_10yr	11.3%	Calculated from inputs
fasting_glucose	9.1%	Median imputation
BMI	1.2%	Median imputation
bp_meds	0.5%	Mode (0)

## Feature Engineering Pipeline

34 engineered features including: - **Derived:** pulse\_pressure, mean\_arterial\_pressure, metabolic\_syndrome\_score - **\*\*Log Transforms:\*\*** log\_cholesterol, log\_glucose, log\_bmi - **Interactions:** age×systolic\_bp, bmi×glucose, age×smoking - **\*\*Categorical:\*\*** age\_group (5 bins), bmi\_category (4 bins) - **\*\*Binary Flags:\*\*** hypertension, high\_cholesterol, high\_glucose, obesity



### 3. Model Architecture

#### Production Model: MLP Classifier



#### Model Specifications

Component	Configuration
Architecture	Multi-Layer Perceptron


Component	Configuration
Hidden Layers	(100, 50) neurons
Activation	ReLU
Output	Softmax (3 classes)
Optimizer	Adam
Regularization	L2 (alpha=0.0001)
Early Stopping	Yes (patience=10)

## 4. Classification Models Comparison

Complete Model Comparison Table

Rank	Model	Type	Accuracy	Precision	Recall	F1-Score	ROC-AUC
1	mlp_binary	MLP	99.61%	99.61%	99.61%	99.61%	0.999
2	final_classifier	MLP	99.25%	99.20%	99.25%	99.22%	0.998
3	stacking_tree_ensemble	Ensemble	99.21%	99.30%	98.94%	99.11%	0.997
4	stacking_lr_ensemble	Ensemble	99.12%	98.99%	98.88%	98.94%	0.996
5	hgb_multiclass_calibrated	HGB	99.08%	99.20%	98.89%	99.04%	0.995
6	mlp_3class	MLP	99.04%	98.79%	98.88%	98.84%	0.994
7	best_real_outcome_model	MLP Deep	98.94%	97.76%	96.11%	96.49%	0.978
8	mlp_multiclass_calibrated	MLP Cal	98.73%	98.67%	98.08%	98.37%	0.991
9	voting_ensemble	Voting	98.60%	98.52%	98.36%	98.44%	0.989
10	svm_binary	SVM	97.81%	97.83%	97.80%	97.81%	0.995
11	rf_binary	RandomForest	97.19%	97.20%	97.19%	97.19%	0.997
12	svm_3class	SVM	96.05%	95.94%	94.75%	95.32%	0.981
13	rf_multiclass_calibrated	RF Cal	95.13%	94.65%	95.22%	94.93%	0.978

Rank	Model	Type	Accuracy	Precision	Recall	F1-Score	ROC-AUC
14	lr_binary	LogReg	94.91%	94.91%	94.92%	94.91%	0.989
15	rf_3class	RandomForest	94.34%	93.45%	94.55%	93.98%	0.975
16	lr_3class	LogReg	91.89%	91.60%	90.30%	90.89%	0.962



## Model Selection Justification

**Selected Production Model:** `final_classifier.pkl` (MLP)


Reasons: 1. **High Accuracy (99.25%)** exceeds the 85% target by significant margin  
2. **Balanced Performance** across all three risk classes  
3. **Fast Inference** (~500ms per prediction)  
4. **Robust to noise** in OCR-extracted data  
5. **Clinical Override Compatible** for edge case handling

---

## 5. Regression Models Comparison

### Regression Model Performance

Model	Type	MAE	RMSE	R <sup>2</sup> Score	Binned Accuracy
hgb_regressor	HistGradientBoosting	<b>0.0075</b>	0.0111	<b>0.992</b>	95.84%
rf_regressor	RandomForest	0.0064	0.0121	0.990	<b>96.45%</b>
risk_regressor_v2	RandomForest	0.0064	0.0121	0.990	96.45%
mlp_regressor	MLP	0.0082	0.0149	0.986	95.35%



### Regression Use Case

Regression models predict **continuous 10-year risk percentage (0-100%)**, which is then binned into risk categories: - LOW: <10% - MODERATE: 10-25% - HIGH: ≥25%

## 6. Hyperparameter Tuning

### Tuning Methodology

- **Method:** RandomizedSearchCV
- **Iterations:** 100 parameter combinations
- **Cross-Validation:** 5-fold stratified

- **Scoring Metric:** F1-score (macro average)
- **Selection Criterion:** Best validation performance

## MLP Classifier Tuned Parameters

Parameter	Search Range	Final Value
hidden_layer_sizes	[(50,), (100,), (100,50), (100,50,25)]	(100, 50)
activation	[relu, tanh]	relu
solver	[adam, sgd]	adam
alpha	[0.0001, 0.001, 0.01]	0.0001
learning_rate	[constant, adaptive]	adaptive
max_iter	[200, 500, 1000]	500
early_stopping	[True, False]	True

## Random Forest Tuned Parameters

Parameter	Search Range	Final Value
n_estimators	[100, 200, 500]	200
max_depth	[10, 20, None]	20
min_samples_split	[2, 5, 10]	5
min_samples_leaf	[1, 2, 4]	2
class_weight	[balanced, None]	balanced

## SVM Tuned Parameters

Parameter	Search Range	Final Value
C	[0.1, 1, 10, 100]	10
kernel	[rbf, linear, poly]	rbf
gamma	[scale, auto]	scale
class_weight	[balanced, None]	balanced



---

## 7. Model Evaluation Metrics

---

### Final Classifier Performance (Test Set)

Metric	LOW	MODERATE	HIGH	Macro Avg
Precision	99.5%	98.8%	99.3%	99.2%
Recall	99.8%	98.5%	99.1%	99.1%
F1-Score	99.6%	98.6%	99.2%	99.2%
Support	1,245	782	392	2,419

### Confusion Matrix Analysis

- **True Positives:** Correctly identified high-risk patients
- **False Negatives:** Minimal - critical for medical safety
- **Specificity:** 99.4% - low false alarm rate

### ROC-AUC Analysis



- **Overall ROC-AUC:** 0.998
- **One-vs-Rest Performance:**
  - LOW vs Rest: 0.999
  - MODERATE vs Rest: 0.996
  - HIGH vs Rest: 0.999


---

## 8. Risk Categorization System

---

### Risk Categories Based on Framingham Score

Category	10-Year Risk	Clinical Interpretation	Action
 <b>LOW</b>	<10%	Low probability of CV event	Maintain healthy lifestyle
 <b>MODERATE</b>	10-25%	Elevated risk, monitoring needed	Lifestyle changes + monitoring

Category	10-Year Risk	Clinical Interpretation	Action
 <b>HIGH</b>	≥25%	Significant risk of CV event	Medical intervention required

## Key Risk Factors (Framingham)

1. **Age** - Risk doubles every 10 years after 45
2. **Blood Pressure** - ≥140/90 mmHg = hypertension
3. **Total Cholesterol** - ≥240 mg/dL = high risk
4. **HDL Cholesterol** - <40 mg/dL = increased risk
5. **Smoking** - Increases risk 2-4x
6. **Diabetes** - CHD risk 2-4x higher

## Risk Class Distribution (Training Data)

Class	Count	Percentage
LOW	7,234	64.1%
MODERATE	2,891	25.6%
HIGH	1,161	10.3%

## 9. Clinical Override Rules

### Problem Addressed

The ML model, trained on Framingham data, may miss young patients with multiple risk factors because: - Young patients (<40): Only 4.1% CHD rate in training data - Model learned "young = low risk" pattern - Clinically dangerous for patients like 32-year-old diabetic smokers

### Clinical Override Implementation

Three safety rules added to `production_model.py` :

### Rule 1: Diabetes Override

```
IF diabetes == 1 AND model_prediction == "LOW":  
    → Override to "MODERATE"
```

**Justification:** Diabetics have 36.7% CHD rate in data

### Rule 2: Young High Metabolic Risk

```
IF age < 50 AND metabolic_score >= 3 AND model_prediction == "LOW":  
    → Override to "MODERATE"
```


**Justification:** Young patients with 3+ risk factors have 15.2% CHD rate

### Rule 3: Extreme Values Safety Net

```
IF systolic_bp >= 180 OR fasting_glucose >= 200:  
    → Override to minimum "MODERATE"
```

**Justification:** Medical emergency values require attention

# Override Impact Analysis



Metric	Value
Total Patients	4,238
Patients Overridden	106 (2.5%)
By Diabetes Rule	7
By Young High Risk	86
By Extreme Values	13
LOW CHD Rate (before)	8.4%
LOW CHD Rate (after)	7.9%

## 10. Saved Models & Pipelines

### Production Model Files

File	Description	Size
models/final_classifier.pkl	Production MLP classifier	568 KB
models/final_classifier_meta.json	Model metadata & feature names	1.2 KB

### Classification Models ( Milestone\_2/models/classification/ )

Model File	Type	Accuracy
final_classifier.pkl	MLP	99.25%
mlp_binary.pkl	MLP	99.61%
mlp_3class.pkl	MLP	99.04%
rf_binary.pkl	Random Forest	97.19%
rf_3class.pkl	Random Forest	94.34%
svm_binary.pkl	SVM	97.81%
svm_3class.pkl	SVM	96.05%
lr_binary.pkl	Logistic Regression	94.91%
lr_3class.pkl	Logistic Regression	91.89%
voting_ensemble.pkl	Voting Ensemble	98.60%
stacking_tree_ensemble.pkl	Stacking	99.21%
stacking_lr_ensemble.pkl	Stacking	99.12%
hgb_multiclass_calibrated.pkl	HGB Calibrated	99.08%
best_real_outcome_model.pkl	MLP (Real CHD)	98.94%

## Regression Models ( Milestone\_2/models/regression/ )







Model File	Type	R <sup>2</sup> Score
hgb_regressor.pkl	HistGradientBoosting	0.992
rf_regressor.pkl	Random Forest	0.990
mlp_regressor.pkl	MLP	0.986

## Preprocessing Pipelines

Component	File	Purpose
OCR Engine	src/production_ocr.py	Extract data from medical images
Feature Engineering	src/production_model.py	Build 34 features from raw data
Clinical Override	src/production_model.py	Apply safety rules
Full Pipeline	src/production_pipeline.py	End-to-end prediction

# 11. Testing & Validation

## Test Cases Validated

Test Report	Expected	Actual	Status
SYN-003 (66yo, multiple risks)	HIGH	HIGH 	Pass
SYN-005 (32yo, healthy)	LOW	LOW 	Pass
SYN-006 (38yo, healthy)	LOW	LOW 	Pass
SYN-007 (58yo, borderline)	MODERATE	MODERATE 	Pass
SYN-009 (78yo, elderly healthy)	MODERATE	MODERATE 	Pass
SYN-010 (32yo, diabetic+smoker)	MODERATE*	MODERATE 	Pass

\*SYN-010 required clinical override (would have been LOW without it)

## Edge Cases Tested






1. **Young patient with severe risks** → Clinical override catches it
2. **Elderly with good vitals** → Correctly MODERATE (age factor)
3. **Missing data fields** → Graceful fallback with defaults
4. **Irrelevant documents (CBC)** → Proper error handling
5. **Poor quality images** → OCR robustness verified

## Error Handling

Scenario	Response
Missing required fields	"Necessary data missing: [field list]"
Invalid document type	"Unable to extract cardiovascular data"
Model loading failure	Fallback to rule-based assessment

# 12. Conclusion

## Achievements

1.  **Exceeded accuracy target** (99.25% vs 85% target)
2.  **Comprehensive model comparison** (18 classification, 4 regression)
3.  **Clinical safety net** implemented with 3 override rules
4.  **Production-ready pipeline** with OCR integration
5.  **Robust error handling** for all edge cases

## Key Technical Contributions

- Advanced feature engineering (34 features from 14 inputs)
- Ensemble methods exploration (voting, stacking)
- Calibrated probability outputs
- Clinical rule integration with ML predictions

## Future Improvements

1. Expand training data with more young patient CHD cases

2. Add explainability features (SHAP values)
  3. Implement confidence intervals for predictions
  4. Mobile app integration via API
- 

**Report Generated:** December 2025

**CardioDetect v2.0** - AI-Powered Cardiovascular Risk Prediction